

Aaron 950™

High Frequency Desiccator with Cut



Distributed By:

Surgisol

Surgical Solutions

201 Trafalgar House, Grenville Place, London NW7 3SA
Tel: 020 8441 650 www.surgisol.com info@surgisol.com

Bovie
AARON
MEDICAL

Aaron 950™ High Frequency Desiccator with Cut

The Aaron 950™, by Bovie®, is the first High Frequency Desiccator with Cut.

It features Cut (60 watts), Blend (60 watts), Coagulation (60 watts), Fulguration / Desiccation (35 watts), and Bipolar (30 watts) modes with six presets for the doctor's favourite settings.

The easy-touch colour-coded membrane pads, on the front panel, allow for quick selection of each setting.

Power output may be increased or decreased via the 3-button handpiece or the rotary power-control.

For maximum flexibility a complete line of disposable and reusable electrodes (blade, ball, needle & loop) and bipolar forceps are available.

Cat # A950-220

DIMENSIONS

Height: 17.8 cm (7")

Width: 26 cm (10.3")

Depth: 11.4 cm (4.5")

Weight: < 3.2 kg (< 7 lbs.)

OUTPUT POWER

Cut: 60 watts

Blend: 60 watts

Coagulation: 60 watts

Fulguration: 35 watts

Bipolar: 30 watts

OUTPUT FREQUENCY

350 – 600 kHz

LINE FREQUENCY

50 – 60 Hz

WARRANTY

one year

STANDARD EQUIPMENT

Comprises: 1 A950 unit, power cord, handpiece (A901), ten nonsterile disposable handpiece drapes (A910), five sterile disposable electrodes (three ES01 blades, one ES02 needle, one ES20 ball), 50 sharp (A804) and 50 blunt (A806) non-sterile disposable dermal tips, reusable grounding cord (A1202C), five disposable grounding pads (A1202), wall mount kit (A837), user's guide



Cat #: S_A950